





Clinician's guide to sample collection for GENITAL sites

Order Box selected on Requisition	APTIMA Unisex Swab 	APTIMA Multitest Swab 	Sterile Orange Top 	Copan Liquid Amies Swab 
<i>Trichomonas (NAAT)</i>	Cervical (F) Urethral (M)	Vaginal Self Collected (F)	Urine (F & M)	
<i>Vaginitis - Initial</i>				Vaginal (F)
<i>Vaginitis – Chronic/recurrent</i>	Cervical (F)	Vaginal Self Collected (F)	Urine (F)	Vaginal (F)
<i>Genital Chlamydia and Gonorrhea (NAAT)</i>	Cervical (F) Urethral (M)	Vaginal Self Collected (F)	Urine (F & M)	
<i>Genital Chlamydia, Gonorrhea and Trichomonas (NAAT)</i>	Cervical (F) Urethral (M)	Vaginal Self Collected (F)	Urine (F & M)	

F = female M= male

F = female

M= male

F = female

M= male